



Bone Instruments



RS-8300

FRIEDMAN
slight curve
very delicate
jaw width 2mm
length 5 1/2"



Micro FRIEDMAN
very delicate
jaw width 1.3mm
length 5 1/2"
straight
curved
very delicate
jaw width 1.0mm
slight curve


RS-8302

RS-8303


RS-8305

RS-8302 Shown


1-800-424-2984
www.roboz.com



RS-8310 Mini Bone Rongeur
jaw width 2mm
jaw length 8mm
overall length 5"




RS-8321 LEMPERT
jaw width 2.5mm
length 6 1/2"
slight curve



RS-8320 Shown

1-800-424-2984
www.roboz.com



RS-8340

Micro BLUMENTHAL
curved
jaw width 2mm
length 4 1/2"



RS-8342


BLUMENTHAL
jaw width 2mm
slight angle
length 6"



RS-8342 Shown

RS-8360

HARTMAN
strong curve
length 5 3/4"



RS-8370

BANE
curved
length 6.5"




RS-8374

ADSON
curved
jaw width 7mm
jaw length 13mm
overall length 8"



75% of Actual






RS-8431

RUSKIN
double action, curved
jaw width 3mm
length 6"

1-800-424-2984
www.roboz.com




RS-8470 Bone Cutter
convex jaws, delicate
length 5"




RS-8480 LITTAUER
bone cutting forceps
straight
length 4 1/2"
RS-8482 length 5"
RS-8484 length 6"



RS-8480 Shown






RS-8490

DEAN
bone cutting forceps
side-cutting
length 6"

RS-8490 Shown




RS-8520

RUSKIN
double action, straight
length 5 1/2"




- RS-8540** LISTON straight, length 5¹/₂"
- RS-8544** straight, length 8¹/₂"
- RS-8555** angled on flat, length 5¹/₂"
- RS-8559** angled on flat, length 8³/₄"



RS-8555 Shown

1-800-424-2984
www.roboz.com



RS-8820

FREER
double ended
sharp and blunt blades
width 5mm
length 7"




RS-8830

MOLT
double ended
sharp and blunt blades
length 7¹/₄"



1-800-424-2984
www.roboz.com



RS-8910
FREER
width 4mm
length 6 1/2"
straight



RS-8910 Shown


RS-9782
VIRCHOW
skull breaker
blade 1/2"
length 3"



RS-9040
RS-9052
SPRATT
length 6 1/2"
oval, size 5/0
oval, size 2



RS-9040 Shown



RS-9110


GLUCK
rib shears
length 8¹/₂"

75% of Actual Size


RS-9120

LANGENBECK
saw
blade length 4¹/₂"
overall length 9"

50% of Actual Size




RS-6280C Round
sizes 1/2 thru 8
hand piece



Please specify size when ordering

RS-6282C Inverted Cone
sizes 33.5 through 39
hand piece




Please specify size when ordering

RS-6306 Trephine
micro dissecting
hand piece
carbide
diameter 1.8mm
RS-6308 diameter 2.3mm
RS-6310 diameter 2.7mm




RS-6308 Shown

1-800-424-2984
www.roboz.com



MICHELE
trepine
graduations 1/4" to 1 1/4"
Shaft length 6 1/4"
Shaft O.D. 3/16"

RS-9200



RS-9200 Shown
80% of Actual Size

1-800-424-2984
www.roboz.com

