

RS-5068 McPHERSON
straight
smooth tying platform
tip 0.3 x 5mm
length 4"



RS-5069 McPHERSON
angled
smooth tying platform
tip 0.4 x 5mm
length 4"



RS-5070 Micro Dissecting
straight
fine tips
length 3 1/4"



RS-5075 Micro Dissecting
cross action
fine tips
length 3 1/2"



RS-5079 Micro Dissecting
smooth jaws, angular
extra fine tips
length 3 1/2"



RS-5090 Micro Dissecting
straight
fine sharp tips
length 4"



RS-5095 Micro Dissecting
angled
fine sharp tips
length 4"



RS-5096 HARMS
micro suturing
straight w/ tying platform
tip width 0.6mm
length 4"





RS-5097 HARMIS
micro suturing
curved w/ tying platform
tip width 0.6mm
length 4"



RS-5100 FOERSTER
straight
serrated tips
tip width 0.5mm
length 4"



RS-5101 FOERSTER
curved
serrated tips
tip width 0.5mm
length 4"



RS-5102 FOERSTER
straight
1x2 teeth
tip width 0.5mm
length 4"




RS-5110 GRAEFE
serrated, straight
tip width 0.6mm
length 2³/₄"



RS-5111 GRAEFE
serrated, curved
tip width 0.6mm
length 2³/₄"



RS-5112 GRAEFE
1x2 teeth, straight
tip width 0.6mm
length 2³/₄"



RS-5113 GRAEFE
1x2 teeth, curved
tip width 0.6mm
length 2³/₄"

RS-5130 Micro Dissecting
serrated, straight
tip width 0.8mm
length 4"

RS-5132 Micro Dissecting
serrated, straight
extra delicate
tip width 0.5mm
length 4"

RS-5135 Micro Dissecting
serrated, slight curve
tip width 0.8mm
length 4"


RS-5136 Micro Dissecting
serrated, slight curve
extra delicate
tip width 0.5mm
length 4"

RS-5137 Micro Dissecting
serrated, full curve
tip width 0.8mm
length 4"

RS-5138 Micro Dissecting
serrated, full curve
extra delicate
tip width 0.5mm
length 4"

RS-5139 Micro Dissecting
serrated, angular
tip width 0.8mm
length 4"





RS-5150 Micro Dissecting
1x2 teeth, straight
tip width 0.8mm
length 4"



RS-5152 Micro Dissecting
1x2 teeth, straight
extra delicate
tip width 0.5mm
length 4"



RS-5153 Micro Dissecting
1x2 teeth, slight curve
tip width 0.8mm
length 4"



RS-5155 Micro Dissecting
1x2 teeth, slight curve
extra delicate
tip width 0.5mm
length 4"



RS-5157 Micro Dissecting
1x2 teeth, full curve
tip width 0.8mm
length 4"



RS-5158 Micro Dissecting
1x2 teeth, angular
tip width 0.8mm
length 4"



RS-5163 ARRUGA
micro embryonic
capsule, curved
length 4"

Please note: Handle in photo slightly different than actual instrument



RS-5165 GILL
slight curve
delicate cross serrations
tip width 0.6mm
length 3"



RS-5166 GILL
slight curve, 1x2 teeth
tip width 0.6mm
length 3"



RS-5168 BOTVIN
serrated tips
tip width 0.5mm
length 3"



RS-5170 BONN
tying platform
1x2 teeth
0.12mm teeth
length 2³/₄"



RS-5172 BONN
tying platform
1x2 teeth
0.12mm teeth
length 3³/₄"



RS-5174 McPHERSON
straight
with tying platform
1 x 2 teeth
0.2mm teeth
length 4"



RS-5175 McPHERSON
angled
with tying platform
1 x 2 teeth
0.2mm teeth
length 4"












RS-5180 BISHOP-HARMON
cross serrated
tip width 0.5mm
length 3¹/₄"



RS-5182 BISHOP-HARMON
1x2 teeth
0.3mm teeth
length 3¹/₄"



	RS-5190	STEVENS 1 x 2 teeth tip width 1mm length 4"
	RS-5204	McCULLOUGH smooth jaws tip width 1.5mm length 4"
	RS-5206	McCULLOUGH cross serrated jaws tip width 1.5mm length 4"
	RS-5210	BONNACOLTO longitudinal serrations cross serrated at tip tip width 1.2mm length 4"
	RS-5211	BRACKEN curved fine cross serrations length 4"
	RS-5228	NUGENT angular, smooth jaws tip width 1.2mm length 4 ^{1/4} "
	RS-5230	Micro-ADSON serrated, delicate tip width 0.5mm length 4 ^{3/4} "
	RS-5231	BROWN ADSON 7x7 teeth side grasping length 4 ^{3/4} "



RS-5232
 Micro-ADSON
 1 x 2 teeth, delicate
 tip width 0.5mm
 length 4³/₄"



ADSON
 length 4³/₄"
 1 x 2 teeth
 with serrations
 tip width 1mm

RS-5233



RS-5235
 without serrations
 tip width 1mm

RS-5233 Shown

RS-5234
 ADSON
 2 x 3 teeth
 tip width 1.3mm
 length 4³/₄"




RS-5236
 ADSON
 cross serrated
 carbide jaws
 length 4³/₄"




RS-5237
 HUDSON (EWALD)
 serrated
 length 4³/₄"



RS-5238
 HUDSON (EWALD)
 length 4³/₄"
 1 x 2 teeth





RS-5242 Shown

RS-5240

RS-5242

RS-5244

SEMKEN
length 5"
serrated
tip width 1.2mm
1x2 teeth
tip width 1.2mm
2x3 teeth
tip width 2mm



RS-5248 Shown

RS-5246

RS-5248

SEMKEN
tip width 1.6mm
length 6"
serrated
1 x 2 teeth



RS-5255

SINGLEY
organ holding forceps
length 6"




RS-5262 Shown

RS-5260

RS-5262

RHOTON
micro suturing
w/ tying platform
length 7"

tip width 0.6mm
tip width 0.3mm



RS-5264 Shown

RS-5263


RS-5264

RHOTON
micro suturing
w/ tying platform
curved
length 5"
length 7"



GERALD
length 7"

- RS-5270 straight, serrated
- RS-5272 straight, 1 x 2 teeth



RS-5270 Shown

GERALD
length 7"

- RS-5271 curved, serrated
- RS-5273 curved, 1 x 2 teeth



RS-5271 Shown

CUSHING
scraper end
length 7"


- RS-5280 serrated
- RS-5282 1 x 2 teeth



RS-5280 Shown


CUSHING
cross serrated
length 7"
carbide jaws

- RS-5284



RS-5310
RS-5312
RS-5314


POTTS-SMITH
serrated
length 7"
length 8¹/₄"
length 9¹/₂"



RS-5310 Shown

RS-5320

POTTS-SMITH
1 x 2 teeth
serrated
length 7"



RS-5320 Shown

RS-5348

Micro Alligator
serrated jaws
tip width 0.8mm
jaw length 4mm
shaft length 3¹/₄"




RS-5358
RS-5359

Micro Dissecting
angled
tip length 11mm
length 6"
serrated jaws
smooth jaws



RS-5358 Shown




- Thumb Dressing
serrated
- RS-8100 length 4 1/2"
- RS-8102 length 5"
- RS-8104 length 5 1/2"
- RS-8106 length 6"
- RS-8108 length 8"
- RS-8110 length 10"
- RS-8112 length 12"



RS-8100 Shown

- Thumb Dressing
delicate, serrated
- RS-8120 length 4 1/2"
- RS-8122 length 5"
- RS-8124 length 5 1/2"
- RS-8126 length 6"




RS-8120 Shown

- RS-8140 Thumb Dressing
serrated
length 3"



- BROWN
8 x 8 teeth
side grasping
- RS-8242 length 6"
- RS-8243 length 8"
- RS-8244 length 10"




RS-8242 Shown

- Splinter Forceps
serrated
fine points
- RS-8150 length 3 1/2"
- RS-8152 length 4"
- RS-8154 length 4 1/2"



RS-8150 Shown





RS-8160 Shown

- RS-8160
- RS-8162
- RS-8164
- RS-8166

Tissue Forceps
 1 x 2 teeth
 length 4 1/2"
 length 5"
 length 5 1/2"
 length 6"



RS-8180 Shown

- RS-8180
- RS-8182
- RS-8184
- RS-8186


Tissue Forceps
 1 x 2 teeth
 delicate
 length 4 1/2"
 length 5"
 length 5 1/2"
 length 6"



RS-8208 Shown

- RS-8208
- RS-8210
- RS-8212
- RS-8214


Tissue Forceps
 2 x 3 teeth
 length 4 1/2"
 length 5"
 length 5 1/2"
 length 6"



RS-8220 Shown

- RS-8220
- RS-8224

Tissue Forceps
 3 x 4 teeth
 length 5"
 length 6"



RS-8240

STILLE
 4 x 5 teeth
 length 5"



RS-8247

HARMON
 4 x 5 teeth
 jaw width 2.5mm
 length 3"



RS-8248

GRAEFE
 with fine teeth
 jaw width 4.5mm
 length 4 1/4"



RS-8249 GRAEFE
locking
jaw width 4.5mm
fine teeth
length 4 1/4"



RS-8254 MOLONEY
serrated
slight curve
length 4 1/2"



RS-8260 KUEHNE
flat jaw
length 4"
angled

RS-8261 straight



RS-8270 RUSSIAN
length 6"



RS-8270 Shown

RS-8280 WILDE
serrated, angular
length 5"



RS-8290

LUCAE
serrated, bayonet
length 5³/₄"



RS-8298

JANSEN
serrated, bayonet
length 8¹/₄"



80% of Actual Size

RS-7560

DeBAKEY
jaw width 2.5mm
length 6"

RS-7562

length 7³/₄"



RS-7560 Shown

RS-7561

DeBAKEY
jaw width 1.5mm
length 6"

RS-7566

length 8"



RS-7561 Shown

