



Provides excellent steam sterilization without moisture build up.

RT-1400S and RT-1400L Shown

These trays are electropolished to provide the toughest, most corrosion-resistant surface and a bright clean look.

Constructed of heavy-duty, 304 18-8 stainless steel, PMP trays are available in many different styles and sizes to meet every sterilization need. End handles allow for safe, easy and secure removal from autoclave.

**Lids sold separately.**

|                                           |                                                                                                                                                                         |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>RT-1350S</b>                           | Instrument Tray<br>dimensions 10x6 <sup>1/2</sup> x3/4"<br>solid                                                                                                        |
| <b>RT-1351 S</b>                          | Instrument Tray<br>dimensions 13 <sup>1/2</sup> x9 <sup>3/4</sup> x3/4"<br>solid                                                                                        |
| <b>RT-1400S<br/>RT-1400P<br/>RT-1400L</b> | Instrument Tray<br>dimensions 9"x5x2"<br>solid<br>perforated<br>lid<br>dimensions 9x5"                                                                                  |
| <b>RT-1401S<br/>RT-1401P<br/>RT-1401L</b> | Instrument Tray<br>dimensions 10x6 <sup>1/2</sup> x2"<br>solid<br>perforated<br>lid<br>dimensions 10x6 <sup>1/2</sup> "                                                 |
| <b>RT-1402S<br/>RT-1402P<br/>RT-1402L</b> | Instrument Tray<br>dimensions 12 <sup>1/4</sup> x7 <sup>3/4</sup> x2 <sup>1/4</sup> "<br>solid<br>perforated<br>lid<br>dimensions 12 <sup>1/4</sup> x7 <sup>3/4</sup> " |
| <b>RT-1500</b>                            | Flash Pan<br>dimensions 9x5x2"                                                                                                                                          |