

Wound Closures





RS-7981

Eye Needles
3/8 circle, cutting edge
size: 1, 2, 3, or 4
dozen

Please specify size when ordering



RS-7982

Eye Needles
1/2 circle, cutting edge
size: 1, 2, 3, or 4
sold by dozen

Please specify size when ordering



RS-7983

Eye Needles
1/2 circle, taper point
size: 1, 2, 3, or 4
sold by dozen

Please specify size when ordering



RS-7985

Regular Surgeon's
3/8 circle, cutting edge
sold by dozen


Please specify size when ordering

RS-7986 Regular Surgeon's
1/2 circle, cutting edge
sold by dozen



Please specify size when ordering

RS-7987 Fine Intestinal Needle
straight, taperpoint
sizes: 1 1/2"
sold by dozen



NON-STERILE

Cat.#	Size	Material		Length	Needle	Packaged
SUT-15-1	5-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool
SUT-15-2	4-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool
SUT-15-3	3-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool
SUT-15-4	2-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool

STERILE


Cat.#	Size	Material		Length	Needle	Packaged
SUT-1234-13	6-0	Silk, Black, Braided	Non-Absorbable	18"	HE-6	Box/12
SUT-1228-13	6-0	Silk, Black, Braided	Non-Absorbable	18"	HE-7	Box/12
SUT-1073-11	6-0	Silk, Black, Braided	Non-Absorbable	18"	C-1	Box/36
SUT-1073-21	5-0	Silk, Black, Braided	Non-Absorbable	18"	C-1	Box/36
SUT-1074-21	5-0	Silk, Black, Braided	Non-Absorbable	18"	C-13	Box/36
SUT-1091-23	5-0	Silk, Black, Braided	Non-Absorbable	18"	P-12	Box/12
SUT-570-31	4-0	Gut, Chromic	Absorbable	18"	C-12	Box/36
SUT-566-31	4-0	Gut, Chromic	Absorbable	18"	C-13	Box/36
SUT-1066-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-22	Box/36
SUT-1073-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-1	Box/36
SUT-1074-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-13	Box/36
SUT-1076-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-14	Box/36
SUT-1203-33	4-0	Silk, Black, Braided	Non-Absorbable	18"	P-13	Box/12
SUT-1123-31	4-0	Silk, Black, Braided	Non-Absorbable	30"	V-20	Box/36
SUT-566-41	3-0	Gut, Chromic	Absorbable	18"	C-13	Box/36
SUT-1074-41	3-0	Silk, Black, Braided	Non-Absorbable	18"	C-13	Box/36
SUT-1076-41	3-0	Silk, Black, Braided	Non-Absorbable	18"	C-14	Box/36
SUT-1123-41	3-0	Silk, Black, Braided	Non-Absorbable	30"	V-20	Box/36

1-800-424-2984
www.roboz.com



NEEDLE CHART FOR SUTURES WITH NEEDLES


REVERSE CUTTING
3/8 CIRCLE




REVERSE CUTTING
1/2 CIRCLE



TAPER
1/2 CIRCLE





RS-7069

Wire Cutting
serrated
angled with notch
length 5"



RS-5950

SPENCER
stitch, very delicate
length 3 1/2"



RS-7074
RS-7076


LITTAUER
stitch
length 4 1/2"
length 5 1/2"

RS-7074 Shown



RS-6400

KALT
with spring catch
length 5 1/2"



RS-6403 PATTON
 straight without lock, serrated
 length 4^{1/2"}
 serrated




RS-6410 CASTROVIEJO
 straight without catch
 length 5^{3/4"}
RS-6412 straight with catch




RS-6412 Shown

RS-6416 CASTROVIEJO
 extra delicate
 0.6mm wide jaws
 to hold 7-0 to 10-0
 sutures
 length 5^{3/4"}
 straight with catch
RS-6417 curved with catch



RS-6417 Shown





RS-6421 Shown

RS-6421

McPHERSON
length 4"
curved
smooth tapered jaws



RS-6430

CASTROVIEJO
round handle with lock
carbide jaws
length 4 1/2"



RS-6437 Shown

RS-6437


Micro Needle Holder
0.6mm tips
length 5 1/8"
curved jaws with lock



RS-6445 Shown

RS-6445

Micro Needle Holder
0.6mm tips
length 7 1/8"
curved jaws with lock



Intracardiac
extra delicate
with carbide jaws
RS-7800 length 5"
RS-7804 length 7"




Note: tungsten carbide inserts have 4,000 teeth per square inch

RS-7820 BABY DERF
large finger rings
length 4^{1/4}"



1-800-424-2984
www.roboz.com





RS-7822
RS-7824

DERF
length 4³/₄"
standard jaws
carbide jaws


RS-7822 Shown



RS-7830
RS-7832
RS-7842

WEBSTER
smooth jaws, delicate
length 4³/₄"
length 5"
serrated carbide jaws
length 5"

RS-7830 Shown




- HALSEY
- length 5"
- RS-7840** smooth jaws
- RS-7841** serrated jaws
- carbide
- RS-7841L** left handed version




RS-7840 Shown

- RS-7850** HEGAR-BAUMGARTNER
- length 5¹/₄"



RS-7860
RS-7862

CRILE-WOOD
length 6"
length 7"




RS-7860 Shown

RS-7880
RS-7882

OLSEN-HEGAR
extra delicate with
combination scissor
length 4^{3/4}"
serrated jaws
extra delicate
smooth



RS-7880 Shown



OLSEN-HEGAR
 needle holder scissors
 combination cross
 serration tip


- RS-7884** length 5 1/2"
- RS-7886** length 6 1/2"



RS-7884 Shown


OLSEN-HEGAR
 needle holder scissors
 combination
 with carbide jaws
 length 4 3/4" TC

- RS-7892 length 4 3/4" TC
- RS-7894** length 5 1/2"
- RS-7896** length 6 1/2"




RS-7894 Shown

Note: tungsten carbide inserts have 2,600 teeth per square inch



RS-7900 METZENBAUM
length 7"
straight




RS-7910 MAYO-HEGAR
length 5"
RS-7912 length 6"
RS-7912L left handed version
RS-7920 with carbide jaws
length 5"
RS-7922 with carbide jaws
length 6"




RS-7910 Shown


Note: tungsten carbide inserts have 2,600 teeth per square inch




RS-7960
BROWN
convex jaw
with fine serrations
for easy suture tying
length 5¹/₄"



RS-7960 Shown





Reflex Wound Clip Products

- RS-9260** Reflex 9 Wound Clip
Applier for Reflex 9mm
clips only

- RS-9262** Reflex 9 Wound Clip
9mm for use with
RS-9260
box of 100
Stainless Steel

- RS-9265** Reflex 9 Wound Clip
9mm for use with
RS-9260
box of 1,000
Stainless Steel

- RS-9250** Reflex 7 Clip Applier
for Reflex 7mm clips
only


- RS-9258** Reflex 7 Clip 7mm
for use with RS-9250
box of 100
Stainless Steel

- RS-9255** Reflex 7 Clip 7mm
for use with RS-9250
box of 1,000
Stainless Steel

-
- RS-9263** Reflex Wound Clip
removing forceps
for 9mm, 7mm and
other brands
length 4"



1-800-424-2984
www.roboz.com



MICHEL
wound clips
packaged 1,000 pcs.
per box

RS-9270 7.5mm

RS-9272 11mm




RS-9290 HEGENBARTH
clip applying forceps
cross action
length 5"



RS-9292 MICHEL
clip applying forceps
length 5"



RS-9294 MICHEL
clip applying and
removing forceps
length 5"



**COMPARISON OF
MILLIMETERS TO INCHES**

MM	INCHES	MM	INCHES
0.33	.013	6	.236
0.5	.020	6.33	.249
0.67	.026	6.5	.256
1	.039	6.67	.263
1.33	.052	7	.275
1.5	.059	7.33	.288
1.67	.066	7.5	.295
2	.079	7.67	.302
2.33	.092	8	.315
2.5	.098	8.33	.328
2.67	.105	8.5	.335
3	.118	8.67	.341
3.33	.131	9	.354
3.5	.138	9.33	.367
3.67	.144	9.5	.374
4	.157	9.67	.380
4.33	.170	10	.394
4.5	.177	10.33	.407
4.67	.184	10.67	.420
5	.197	11	.433
5.33	.210	11.33	.446
5.5	.216	11.67	.459
5.67	.223	12	.472

**COMPARISON OF
INCHES TO MILLIMETERS**

INCHES	MM	INCHES	MM
1/32	0.79	17/32	13.49
1/16	1.59	9/16	14.29
3/32	2.38	19/32	15.08
1/8	3.18	5/8	15.88
5/32	3.97	21/32	16.67
3/16	4.76	11/16	17.46
7/32	5.56	23/32	18.26
1/4	6.35	3/4	19.05
9/32	7.14	25/32	19.84
5/16	7.94	13/16	20.64
11/32	8.73	27/32	21.43
3/8	9.53	7/8	22.22
13/32	10.32	29/32	23.02
7/16	11.11	15/16	23.81
15/32	11.91	31/32	24.61
1/2	12.7	1"	25.4

**COMPARISON OF
INCHES TO
CENTIMETERS**

INCHES	CM
1	2.54
2	5.08
3	7.62
4	10.16
5	12.7
6	15.24
7	17.78
8	20.32
9	22.86
10	25.40
11	27.94
12	30.48
13	33.02
14	35.56
15	38.10

**COMPARISON OF
CENTIMETERS TO
INCHES**

CM	INCHES
1	.3937
2	.7874
3	1.1811
4	1.5748
5	1.9685
6	2.3622
7	2.7559
8	3.1496
9	3.5433
10	3.937
20	7.8740
30	11.811
40	15.748

