

32 ga micro injection needles uses:

- portal vein injections
- thymus injections
- bile duct injections
- blood/fluid extractions

Advantages:

- potential for clotting is reduced
- reduces trauma
- greater control of injection site
- greater volume control



IN-900


Micro Injection Needle
 .009" OD x .003" ID
 32 ga
 needle length 10mm
 use with IN-910

straight



IN-905

45 deg angle



IN-910


IN-910

Syringe
 plastic
 Volume 1mm
 for use with IN-900
 and IN-905 needles

Photo Not Available For IN-940

IN-940

Micro Injection Needle
 straight
 .009" OD x .003" ID
 32 ga
 needle length 12mm
 use with IN-950



IN-950

IN-950

Syringe
 glass
 5 microliters
 for use with IN-940
 needles

