

**RS-6470**      Micro Vascular Clip  
 vessel diameter 1mm  
 occlusion pressure 30-35gr  
 length 8mm



**RS-6472**      Micro Vascular Clip  
 vessel diameter 1-2mm  
 occlusion pressure 35-40gr  
 length 11mm



**RS-6474**      Micro Vascular Clip  
 vessel diameter 2-3mm  
 occlusion pressure 50-60gr  
 length 14mm



**RS-6488**      Micro Vascular Approximator  
 with moveable clips  
 vessel diameter 1mm  
 occlusion pressure 30-35gr  
 length 8mm



**RS-6490**      Micro Vascular Approximator  
 with moveable clips  
 vessel diameter 1-2mm  
 occlusion pressure 35-40gr  
 length 11mm



**RS-6492**      Micro Vascular Approximator  
 with moveable clips  
 vessel diameter 2-3mm  
 occlusion pressure 50-60gr  
 length 14mm



**RS-6496**      Micro Clip Setting  
 Forceps  
 length 4<sup>3</sup>/<sub>4</sub>"



To avoid damage and for ease of use, please use the micro vascular clips with these clip setting forceps.

