

New Disposable Clippers from Roboz Surgical

- Removes all types of head and body hair, wet or dry.
- Unique design pivoting head rotates easily to allow flexibility of clipping with either a pushing or a pulling (raking) motion.
- The cutting blade is on top of the stationary blade, away from the skin, so there's minimal risk of nicks or cuts.
- Single-use blades help reduce the risk of cross-contamination and are easy to snap on for use and pop off for disposal.
- Rechargeable clipper provides approximately 60 minutes of run time, from a single 8-hour recharge.
- Unique design is easy to clean and disinfect, wipe with a disinfectant or rinse under running water.

Cordless Clipper and Charger

RC-5902

\$208.00

Pack of 5 Disposable Blades

RC-5902RB

\$22.20

Call **1-800-424-2984** for Details