

Clips and Clamps



		RS-5420	Micro Clip straight 70gr pressure clip width 0.75mm jaw length 5mm
		RS-5422	Micro Clip straight 70gr pressure clip width 0.75mm jaw length 6mm
		RS-5424	Micro Clip straight 70gr pressure clip width 1mm jaw length 8mm
		RS-5426	Micro Clip straight 70gr pressure clip width 1.5mm jaw length 10mm
		RS-5431	Micro Clip slight curve 70gr pressure clip width 1mm jaw length 6mm
		RS-5433	Micro Clip strong curve 70gr pressure clip width 1mm jaw length 6mm
		RS-5435	Micro Clip 45° angle 70gr pressure clip width 1mm jaw length 6mm


Please use clip applying forceps with these micro clips to avoid damage and for ease of use.



RS-5410
Clip Applying Forceps
for micro clips
length 5¹/₂"




RS-5440
Micro Clip
applying forceps
with lock
length 5³/₄"



RS-5450

SCHWARTZ clip applying forceps for temporary clips length 8"



75% of actual size

RS-5452

SCHWARTZ temporary clip straight, smooth 795gr pressure jaw width 1.7mm jaw length 8mm clip length 1"



RS-5457

SCHWARTZ temporary clip slight angle, smooth 795gr pressure jaw width 1.7mm jaw length 8mm clip length 1"



RS-5459

SCHWARTZ temporary clip strong angle, smooth 795gr pressure jaw width 1.7mm jaw length 8mm clip length 1"



RS-6470 Micro Vascular Clip
 vessel diameter 1mm
 occlusion pressure 30-35gr
 length 8mm



RS-6472 Micro Vascular Clip
 vessel diameter 1-2mm
 occlusion pressure 35-40gr
 length 11mm



RS-6474 Micro Vascular Clip
 vessel diameter 2-3mm
 occlusion pressure 50-60gr
 length 14mm



RS-6488 Micro Vascular Approximator
 with moveable clips
 vessel diameter 1mm
 occlusion pressure 30-35gr
 length 8mm




RS-6490 Micro Vascular Approximator
 with moveable clips
 vessel diameter 1-2mm
 occlusion pressure 35-40gr
 length 11mm




RS-6492 Micro Vascular Approximator
 with moveable clips
 vessel diameter 2-3mm
 occlusion pressure 50-60gr
 length 14mm



RS-6496 Micro Clip Setting
 Forceps
 length 4³/₄"



To avoid damage and for ease of use, please use the micro vascular clips with these clip setting forceps.





RS-5470

DIEFFENBACH
 serrefine
 straight, very delicate
 300gr pressure
 length 1 1/8"



RS-5471

DIEFFENBACH
 serrefine
 curved, very delicate
 300gr pressure
 length 1 1/8"



RS-7438

DIETRICH
 micro bulldog clamp
 straight
 80gr pressure
 jaw length 8mm
 jaw width 1.2mm
 length 2"



RS-7439

DIETRICH
 micro bulldog clamp
 curved
 80gr pressure
 jaw length 8mm
 jaw width 1.2mm
 length 2"



RS-7440 JOHNS HOPKINS
 bulldog clamp
 straight
 700gr pressure
 length 1 1/2"



RS-7441 JOHNS HOPKINS
 bulldog clamp
 curved
 700gr pressure
 length 1 1/2"



RS-7442 JOHNS HOPKINS
 bulldog clamp
 straight
 800gr pressure
 length 2"



RS-7443 JOHNS HOPKINS
 bulldog clamp
 curved
 800gr pressure
 length 2"



RS-7444 JOHNS HOPKINS
 bulldog clamp
 straight
 800gr pressure
 length 2 1/2"




RS-7445 JOHNS HOPKINS
 bulldog clamp
 curved
 800gr pressure
 length 2 1/2"



RS-7446 JOHNS HOPKINS
 bulldog clamp
 straight
 1000gr pressure
 length 3"



RS-7447 JOHNS HOPKINS
 bulldog clamp
 curved
 1000gr pressure
 length 3"





RS-7500

DeBAKEY
atraumatic bulldog clamp
straight
400-500gr pressure
jaw length 2cm
overall length 7.5cm



RS-7511

DeBAKEY
atraumatic bulldog clamp
curved
400-500gr pressure
jaw length 2cm
overall length 7.5cm



RS-7530 GLOVER
atraumatic bulldog clamp
straight
tension adjusting screw
240-260gr pressure
jaw length 2cm
overall length 7cm



RS-7541 GLOVER
atraumatic bulldog clamp
curved
with tension adjusting screw
240-260gr pressure
jaw length 2cm
overall length 7cm





RS-7760

JONES
towel clamp
cross action
light weight
length 2¹/₄"



RS-7784 Shown

RS-7780

BACKHAUS
towel clamp
length 3¹/₂"

RS-7784

length 5¹/₄"

