

RS-5420



Micro Clip
straight
60gr pressure
clip width 0.8mm
jaw length 5mm

RS-5422



Micro Clip
straight
85gr pressure
clip width 0.75mm
jaw length 6mm

RS-5424



Micro Clip
straight
100gr pressure
clip width 1mm
jaw length 8mm

RS-5426



Micro Clip
straight
125gr pressure
clip width 1.5mm
jaw length 10mm

RS-5431



Micro Clip
slight curve
85gr pressure
clip width 1mm
jaw length 6mm

RS-5433



Micro Clip
strong curve
85gr pressure
clip width 1mm
jaw length 6mm

RS-5435



Micro Clip
45° angle
85gr pressure
clip width 1mm
jaw length 6mm

