

Micro Dissecting Needles And Insect Pins



- **Adjustable length**
- **Easy to handle**
- **Durable material**
- **Finer and less costly than tweezers**

Sharp Point Tungsten Needles

These fine needles made of hard tungsten have a 1 μ m tip diameter. They can be used for micro dissection at the finest level including the separation of cells.

The needles are 2" long for easy handling. Insert the needles deeper into the handle to gain lateral stability or leave a longer piece out to increase flexibility.

Strong Point Micro Dissecting Needles







These micro dissecting needles have super fine yet strong tips with 10 μ m diameter. They are suited for dissections at the microscopic scale which require increased lateral strength. Strong point needles are available with straight and angled tips.

Insect Pins And Minutiens

Two kinds of fine, hand made pins are available. Stainless steel pins with nylon heads are typically used in humid environments. Black enameled pins offer extra protection against corrosion and also reduce glare. For the mounting of very small insects or other specimen, 0.25-0.35 mm diameter insect pins can be used as well as headless minutiens.



Available Products

Product	Item #	Quantity	Rod diameter mm	Length	Price \$	
Micro Dissecting Needles						
Micro Dissecting Needle Holders 	RS-6060	1		5 1/4"	24.00	
	RS-6061	1		4 3/4"	17.00	
	RS-6062	1		4 in	12.00	
Sharp Point Tungsten Needles 	RS-6063	10	0.125	2"	87.00	
	RS-6064	10	0.25	2"	87.00	
	RS-6065	10	0.5	1.5"	87.00	
Strong Point Micro Dissection Needles 	Angled 45 degrees	RS-6066	5	0.5	2"	105.00
	Angled 90 degrees	RS-6067	5	0.5	2"	140.00
	Straight	RS-6068	5	0.5	2"	140.00
	Angled 45 degrees	RS-6066-1	1	0.5	2"	35.00
	Angled 45 degrees	RS-6067-1	1	0.5	2"	46.70
	Angled 90 degrees	RS-6068-1	1	0.5	2"	46.70
Insect Pins And Minutien Pins						
Black Enameled Insect Pins 	RS-6081-25	100	0.25	38 mm	8.00	
	RS-6081-30	100	0.3	38 mm	8.00	
	RS-6081-35	100	0.35	38 mm	8.00	
	RS-6081-40	100	0.4	38 mm	8.00	
	RS-6081-45	100	0.45	38 mm	8.00	
	RS-6081-50	100	0.5	38 mm	8.00	
	RS-6081-55	100	0.55	38 mm	8.00	
	RS-6081-60	100	0.6	38 mm	8.00	
	RS-6081-65	100	0.65	38 mm	8.00	
	RS-6081-70	100	0.7	52 mm	8.00	
Stainless Steel Insect Pins 	RS-6082-25	100	0.25	38 mm	9.00	
	RS-6082-30	100	0.3	38 mm	9.00	
	RS-6082-35	100	0.35	38 mm	9.00	
	RS-6082-40	100	0.4	38 mm	9.00	
	RS-6082-45	100	0.45	38 mm	9.00	
	RS-6082-50	100	0.5	38 mm	9.00	
	RS-6082-55	100	0.55	38 mm	9.00	
	RS-6082-60	100	0.6	38 mm	9.00	
	RS-6082-65	100	0.65	38 mm	9.00	
	RS-6082-70	100	0.7	52 mm	9.00	
Minutien pins 	stainless steel	RS-6083-10	500	0.1	12 mm	26.00
	stainless steel	RS-6083-15	500	0.15	12 mm	26.00
	stainless steel	RS-6083-20	500	0.2	12 mm	26.00
	black enameled	RS-6084-15	500	0.15	12 mm	28.00

