



RS-5620

VANNAS  
straight  
cutting edge 3mm  
comb. tip width 0.2mm  
overall length 3"



RS-5621

VANNAS  
curved  
cutting edge 3mm  
comb. tip width 0.2mm  
overall length 3"



RS-5650

CASTROVIEJO  
straight, sharp  
cutting edge 9mm  
comb. tip width 0.3mm  
overall length 3<sup>1</sup>/<sub>2</sub>"



RS-5658

CASTROVIEJO  
angular, sharp  
cutting edge 9mm  
comb. tip width 0.3mm  
overall length 3<sup>1</sup>/<sub>2</sub>"



RS-5668

CASTROVIEJO  
angular, blunt  
cutting edge 9mm  
comb. tip width 0.3mm  
overall length 3<sup>1</sup>/<sub>2</sub>"



RS-5671

CASTROVIEJO  
curved, sharp  
cutting edge 10mm  
comb. tip width 0.33mm  
overall length 3<sup>3</sup>/<sub>4</sub>"

