

**RS-5600**

McPHERSON-VANNAS
straight, sharp
cutting edge 5mm
comb. tip width 0.3mm
overall length 3"

**RS-5601**

McPHERSON-VANNAS
curved, sharp
cutting edge 5mm
comb. tip width 0.3mm
overall length 3"

**RS-5602**

McPHERSON-VANNAS
straight, sharp
cutting edge 5mm
comb. tip width 0.2mm
overall length 3"

**RS-5603**

McPHERSON-VANNAS
curved, sharp
cutting edge 5mm
comb. tip width 0.2mm
overall length 3"

**RS-5604**

McPHERSON-VANNAS
straight
sharp needle points
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"

**RS-5605**

McPHERSON-VANNAS
curved
sharp needle points
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"



RS-5630 McPHERSON-VANNAS
 straight
 cutting edge 6mm
 comb. tip width 0.3mm
 overall length 4"



RS-5631 McPHERSON-VANNAS
 curved
 cutting edge 6mm
 comb. tip width 0.3mm
 overall length 4"



RS-5610 VANNAS
 straight
 cutting edge 3mm
 comb. tip width 0.3mm
 overall length 3"




RS-5611 VANNAS
 curved
 cutting edge 3mm
 comb. tip width 0.3mm
 overall length 3"



RS-5618 VANNAS
 angled on edge
 cutting edge 3mm
 comb. tip width 0.3mm
 overall length 3"



RS-5619 VANNAS
 angled on flat
 cutting edge 3mm
 comb. tip width 0.3mm
 overall length 3"





RS-5620

VANNAS
straight
cutting edge 3mm
comb. tip width 0.2mm
overall length 3"



RS-5621

VANNAS
curved
cutting edge 3mm
comb. tip width 0.2mm
overall length 3"



RS-5650

CASTROVIEJO
straight, sharp
cutting edge 9mm
comb. tip width 0.3mm
overall length 3¹/₂"



RS-5658

CASTROVIEJO
angular, sharp
cutting edge 9mm
comb. tip width 0.3mm
overall length 3¹/₂"




RS-5668

CASTROVIEJO
angular, blunt
cutting edge 9mm
comb. tip width 0.3mm
overall length 3¹/₂"



RS-5671

CASTROVIEJO
curved, sharp
cutting edge 10mm
comb. tip width 0.33mm
overall length 3³/₄"



RS-5673 CASTROVIEJO
 curved, blunt
 cutting edge 10mm
 comb. tip width 0.33mm
 overall length 3¹/₂"



RS-5674 CASTROVIEJO
 13mm blades
 angled on flat, blunt
 overall length 4¹/₂"



RS-5675 NOYES
 slightly curved, sharp
 cutting edge 13mm
 comb. tip width 0.4mm
 overall length 4¹/₂"




RS-5676 NOYES
 straight, sharp
 cutting edge 13mm
 comb. tip width 0.4mm
 overall length 4¹/₂"



RS-5677 NOYES
 straight, sharp-blunt
 overall length 4¹/₂"



RS-5678 NOYES
 straight, blunt-blunt
 overall length 4¹/₂"



RS-5680

Micro Dissecting
straight, 8mm blades
flat handle
comb. tip width 0.33mm
overall length 6"

**RS-5681**

Micro Dissecting
curved, 8mm blades
flat handle
comb. tip width 0.33mm
overall length 6"

**RS-5692**

Micro Dissecting
straight, 6mm blades
round handle
comb. tip width 0.33mm
overall length 6"



RS-5693 McCLURE
 Micro Dissecting
 curved, 6mm blades
 round handle
 comb. tip width 0.33mm
 overall length 6"




RS-5788 NOYES
 angled on flat
 12mm blades
 length 4 3/4"



RS-5802 DeWECKER
 sharp-sharp
 7mm blades
 length 4 1/4"



1-800-424-2984
www.roboz.com



Ordering Information

Phone To Order by phone call Toll Free **1-800-424-2984**
Our experienced customer support staff will courteously and efficiently assist you. Our business hours are:
8:30 a.m. to 5:30 p.m. Eastern Time

Fax Fax your purchase orders directly to:
Roboz Surgical Customer Service.
1-888-424-3121

Mail Mail your purchase orders to:
Roboz Surgical Instrument Co., Inc.
P.O. Box 10710
Gaithersburg, Maryland 20898-0710

Credit Card Roboz accepts VISA, Mastercard and American Express   

After Hours Our normal business hours are: 8:30 a.m. to 5:30 p.m. Eastern Time
Holidays & Weekends Monday through Friday. Outside normal hours please leave your name on our voice mail and a customer service representative will return your call the next business day.

Shipping & Handling We normally ship UPS two day service for a flat rate fee.
Next day service is available for an additional fee.
Shipping point is Gaithersburg, Maryland.

Returns To return any of our products please call us at 1-800-424-2984. Our customer service staff will arrange for authorization and assist you. Please use care in packaging delicate items. Also note a 10% restocking charge may apply.

Terms & Conditions **Net 30 Days**
A 1.5% charge for overdue accounts may apply.
Delivery 2-5 days ARO
F.O.B. Gaithersburg, MD