

RS-6403 PATTON
 straight without lock, serrated
 length 4^{1/2}"
 serrated




RS-6410 CASTROVIEJO
 length 5^{3/4}"
 straight without catch
RS-6412 straight with catch




RS-6412 Shown

RS-6416 CASTROVIEJO
 extra delicate
 0.6mm wide jaws
 to hold 7-0 to 10-0
 sutures
 length 5^{3/4}"
 straight with catch
RS-6417 curved with catch



RS-6417 Shown





RS-6421 Shown

RS-6421

McPHERSON
length 4"
curved
smooth tapered jaws



RS-6430

CASTROVIEJO
round handle with lock
carbide jaws
length 4 1/2"



RS-6437 Shown

RS-6437


Micro Needle Holder
0.6mm tips
length 5 1/8"
curved jaws with lock



RS-6445 Shown

RS-6445

Micro Needle Holder
0.6mm tips
length 7 1/8"
curved jaws with lock





RS-6400

KALT
with spring catch
length 5^{1/2}"

