



Provides excellent steam sterilization without moisture build up.
RT-1400S and RT-1400L Shown

These trays are electropolished to provide the toughest, most corrosion-resistant surface and a bright clean look.

Constructed of heavy-duty, 304 18-8 stainless steel, PMP trays are available in many different styles and sizes to meet every sterilization need. End handles allow for safe, easy and secure removal from autoclave.

Lids sold separately.

RT-1350S	Instrument Tray dimensions 10x6 ^{1/2} x3/4" solid
RT-1351 S	Instrument Tray dimensions 13 ^{1/2} x9 ^{3/4} x3/4" solid
RT-1400S RT-1400P RT-1400L	Instrument Tray dimensions 9"x5x2" solid perforated lid dimensions 9x5"
RT-1401S RT-1401P RT-1401L	Instrument Tray dimensions 10x6 ^{1/2} x2" solid perforated lid dimensions 10x6 ^{1/2} "
RT-1402S RT-1402P RT-1402L	Instrument Tray dimensions 12 ^{1/4} x7 ^{3/4} x2 ^{1/4} " solid perforated lid dimensions 12 ^{1/4} x7 ^{3/4} "
RT-1500	Flash Pan dimensions 9x5x2"