

RS-8290

LUCAE
serrated, bayonet
length 5³/₄"



RS-8298

JANSEN
serrated, bayonet
length 7³/₄"



RS-7560

RS-7562

DeBAKEY
jaw width 2.5mm
length 6"
length 7³/₄"



RS-7560 Shown

RS-7561


RS-7566

DeBAKEY
jaw width 1.5mm
length 6"
length 8"




RS-7561 Shown



**RS-7341**


SCHNIDT
slight curve
length 7¹/₂"

75% of Actual Size

**RS-7580
RS-7581**

DeBAKEY
bulldog clamps
ring handle
length 4³/₄"
straight
curved

RS-7580 Shown

**RS-7660**

COOLEY-DERRA
anastomosis clamp
5mm calibrations on jaws
length 6¹/₄"
small clamp

RS-7660 Shown

