

**RS-7800** Intracardiac  
extra delicate  
with carbide jaws  
length 5"  
**RS-7804** length 7"




Note: tungsten carbide inserts have 4,000 teeth per square inch

**RS-7820** BABY DERF  
large finger rings  
length 4<sup>1/4</sup>"



**1-800-424-2984**  
[www.roboz.com](http://www.roboz.com)





**RS-7822**  
**RS-7824**

DERF  
length 4<sup>3</sup>/<sub>4</sub>"  
standard jaws  
carbide jaws


RS-7822 Shown



**RS-7830**  
**RS-7832**  
**RS-7842**

WEBSTER  
smooth jaws, delicate  
length 4<sup>3</sup>/<sub>4</sub>"  
length 5"  
serrated carbide jaws  
length 5"

RS-7830 Shown




- HALSEY
- length 5"
- RS-7840** smooth jaws
- RS-7841** serrated jaws
- carbide
- RS-7841L** left handed version




RS-7840 Shown

- RS-7850** HEGAR-BAUMGARTNER
- length 5<sup>1</sup>/<sub>4</sub>"



**RS-7860**  
**RS-7862**

CRILE-WOOD  
length 6"  
length 7"




RS-7860 Shown

**RS-7880**  
**RS-7882**

OLSEN-HEGAR  
extra delicate with  
combination scissor  
length 4<sup>3/4</sup>"  
serrated jaws  
extra delicate  
smooth



RS-7880 Shown



OLSEN-HEGAR  
 needle holder scissors  
 combination cross  
 serration tip


- RS-7884** length 5 1/2"
- RS-7886** length 6 1/2"



RS-7884 Shown


OLSEN-HEGAR  
 needle holder scissors  
 combination  
 with carbide jaws  
 length 4 3/4" TC

- RS-7892 length 4 3/4" TC
- RS-7894** length 5 1/2"
- RS-7896** length 6 1/2"




RS-7894 Shown

Note: tungsten carbide inserts have 2,600 teeth per square inch



**RS-7900** METZENBAUM  
length 7"  
straight




**RS-7910** MAYO-HEGAR  
length 5"  
**RS-7912** length 6"  
**RS-7912L** left handed version  
**RS-7920** with carbide jaws  
length 5"  
**RS-7922** with carbide jaws  
length 6"




RS-7910 Shown

Note: tungsten carbide inserts have 2,600 teeth per square inch



**RS-7960**  
BROWN  
convex jaw  
with fine serrations  
for easy suture tying  
length 5<sup>1</sup>/<sub>4</sub>"



RS-7960 Shown

